

LIST OF SECONDARY ANTIBODIES AND OTHER CONSUMABLES

Secondary Antibodies:

Donkey-anti-Rabbit Alexa 488---- 20 μ L @ R200 per aliquot
Donkey-anti-Mouse Alexa 488 ---- 20 μ l @ R200 per aliquot
Donkey-anti-Goat Alexa 488 -----20 μ l @ R200 per aliquot
Donkey anti-Rabbit Cy3 ----- 20 μ L @ R200 per aliquot
Donkey anti-Mouse Cy3 ----- 20 μ L @ R200 per aliquot
Donkey anti-Goat CY3 ----- 20 μ l @ R180 per aliquot
Donkey-anti-Mouse Alexa 647 ----- 20 μ l @ R250 per Aliquot
Donkey-anti-Rabbit Alexa 647 ----- 20 μ l @ R220 per Aliquot
Phalloidin_Alexa 488(Green) 50 μ l ----- R360 per aliquot

Organelle Markers:

Yellow Lyso-Tracker ----- For Free
Green Lyso-Tracker ----- 5 μ l @ R32 per Aliquot
----- 50 μ l @ R320 per Aliquot
Green ER-Tracker ----- 5 μ l @ R260 per Aliquot

Others:

- 1) Mowiol mounting medium with anti-fade @ R10 per 1ml
- 2) Hoechst nuclear stain @ R20 per eppe
- 3) Coverslip glass bottom 4 compartment T/C dishes for live cell viewing ----- R180 each

- 4) Coverslip glass bottom single compartment non T/C dishes for live cell viewing ----- R75 each
 - 5) Normal Glass slides and coverslips for microscopy (22x22mm;22x40mm;22x50mm;22x60mm)
 - 6) Superfrosted glass slides @ R7 per slide.
-